

Article Datasheet

Art.No.: 81793

**PH-neutral Rinse Solution, 200ML Bottle
WITH STERILE PHOSPHATE BUFFER SOLUTION (4.9%)**

Article Description:

- in case of cauterization
- 200ml bottle with sterile phosphate buffer solution (4.9%) for neutralizing acids and other substances
- suited with ergonomically formed eye- and dust-protection
- for one-time-use of as refill for eye wash station
- time of use: about 2 minutes
- Packing Unit: 10



Catalogue Page: T 115